


# Classic Folding Wall-Mounted Exam Table Owner's Manual

---



**Model Number:** 12572-00-GZAADG  
**Former Model Number:** 102600-00

- **Two-sided Work Surface**
- **Folds up to conserve space**
- **Easy to use lock/release handle**
- **Stainless Steel Construction - Designed to Last**

#### **SSCI Contact Information**

Contact SSCI Customer Service by mail, telephone or fax from 8:30am to 5:00pm, Central Time, Monday through Friday and closed holidays.

Address: Suburban Surgical Co., Inc. 275 Twelfth street Wheeling, Illinois 60090

Telephone: Illinois-(847)537-9320, ext. 3518 Toll Free-(800)323-7366

Fax: (847)537-9061

Web: [www.suburbansurgical.com](http://www.suburbansurgical.com)

E-mail: [customerservice@suburbansurgical.com](mailto:customerservice@suburbansurgical.com)

Form No. 702710  
Revised 12/05/13

### Care and Cleaning of Stainless Steel Introduction

Stainless steel is steel alloyed with chromium to make it highly resistant to stain, rust, and corrosion. **Note:** This does NOT mean that stainless steel will never rust or corrode. Science has not yet developed a steel which is completely stainless or corrosion PROOF. The type of stainless steel and finish selected by SSCI for this product is the best available for the intended use.

#### Cleaning and Cleansers

The basic rule of thumb is to use the mildest cleaning procedure that will do the job effectively. Always rinse thoroughly with clear water, and dry completely. Frequent cleaning will prolong the service life of stainless steel equipment and will help maintain a bright, pleasing appearance. Ordinary deposits of waste and fluids can usually be removed with soap and water. More stubborn deposits or tightly adhering debris may require harder scrubbing. They also may possibly require the use of commercial cleaning products acceptable for use on metal surfaces. When using any cleaning agent, rub in the direction of the polish lines or "grain" of the metal. For high luster finishes, clean soft cloths or pads should be used. If especially rough cleaning is necessary, use "stainless steel" wool, nylon or plastic scrubbers. Test these scrubbers in an inconspicuous area first to be sure they do not mark or scratch the stainless steel finish. Minor scale build-up and some hard water spotting may be removed by washing with some vinegar, followed by a neutralizing rinse with clear water. A thorough drying with a soft cloth should follow. For heavy deposits of scale, 5% oxalic acid (use warm), 5-15% sulfuric acid, or 5-10% phosphoric acid may be used. Always follow with a neutralizing rinse of clean water and a thorough drying. The carts are fully cart-washable.

#### Deodorizing Agents, Disinfectants & Sanitizers

The large selection of brands and combinations of chemicals available for deodorizing, disinfecting and sanitizing is staggering. Select one or more agents for use in your facility only after weighing in all the benefits claimed by each product. Often this choice is made without adequate consideration of the effects these agents may produce on equipment or furnishings.

Avoid prolonged use of chlorides (such as chlorine bleach), bromides, iodides and thiocyanate on stainless steel surfaces as these chemicals will cause pitting, corrosion and metal discoloration. Allowing salty solutions to evaporate and dry on stainless steel may also contribute to corrosive conditions. In summary, select chemical deodorizers, disinfectants and/or sanitizers only after weighing in all possible outcomes and known adverse effects.

**CAUTION: Before selecting a chemical to employ in your facility, review label statements regarding use with metals (stainless steel). Always consult the chemical supplier if there are any doubts.**

**CAUTION: The warranty for this product is void if the care and cleaning instructions provided in this manual are not followed.**

**Chapter 1 – General Information**

Introduction .....	4
About this Manual .....	4
Information and Safety Notices .....	4
Note: .....	4
CAUTION: .....	4
Parts Ordering Procedure.....	5
Unpacking and Inspection.....	6
Returning the Cart for Repairs .....	6
RMA Numbers .....	6
Packing and Shipment .....	6
Warranty .....	6

**Chapter 2 – Assembly & Installation**

General.....	7
Pre-assembly preparation.....	7
Assembly Procedures.....	7
Tools & Supplies Required.....	7
Drilling the Mounting Holes.....	7
Mounting the Table to the Wall.....	8

**Chapter 3 – Folding & Unfolding the Table**

Folding the Table.....	10
Unfolding the Table.....	10

**Chapter 4 – Accessory (Tie-down Brackets)**

Overview.....	11
Tools Required .....	11
Procedure.....	11

**Chapter 5 – Cleaning**

Introduction.....	12
Cleaning Procedures.....	12

## Chapter 1- General Information

### Introduction

SSCI's Folding Wall-Mounted Exam Table is the ideal solution in confined spaces where an extra work surface is needed. It quickly folds up, out of the way against the wall when not in use. When folded, it extends less than five inches (12.70 cm) from the wall.

The easy to use lock/release handle lets you lower the table when needed and locks the table in either the extended or folded position. The rigid stainless steel frame and support bracket mount the table securely to the wall.

Edges are flanged down and hemmed to help eliminate sharp corners, and a Melamine finished undersurface repels moisture for easy care and long life. The table can be mounted at any convenient height to meet your needs.

### About this Manual

Every attempt has been made to insure that the information in this manual is correct and complete. SSCI, however, always welcomes our customer's suggestions for improvements to our products and associated publications.


### Information & Safety Notices

Throughout this manual you will find text under the headings Note: & CAUTION: The text followed after "Note:" will assist you with additional information about the subject being discussed. The text followed after "CAUTION:" is there to alert you to potentially hazardous conditions which, if ignored or mishandled, could result in injury to yourself or damage to the equipment.

For Example:

**Note:** This does NOT mean that stainless steel will never rust or corrode. Science has not yet developed a steel which is completely stainless or corrosion PROOF.

**CAUTION:** To reduce the possibility of injury, we recommend that unpacking, assembling, installations and replacement operations involving these components be done by at least two people.



**Note:** The Classic Folding Wall Mounted Exam Table (P/N 102600-00) is shipped as a single, pre-assembled unit.

**Note:** Table-to-wall mounting hardware is not included.

**CAUTION:** This table is designed to carry weights up to, but not exceeding, 200-lbs. Placing weights greater than 200-lbs on the table can damage the table or create a hazardous situation.

#### Parts Ordering Procedure

Order new equipment, accessories and replacement parts directly through SSCI Customer Service. You can order by mail, telephone, or fax. Refer to SSCI Contact Information in front for address, Telephone, and fax numbers. When ordering, please provide the following information:

- Your name
- Company name
- Company account number
- Telephone number
- Fax number
- e-mail address
- Shipping address
- Billing address (if different from shipping address)
- Names, part numbers, and quantities of items being ordered
- Credit card number and expiration date, or other payment information
- Preferred method of shipment
- Information on whether the items are required on a normal or urgent basis

## Unpacking and Inspection

**CAUTION: Unpacking the exam table is not difficult. However, it is heavy and we recommend that unpacking be done by at least two people.**

If the shipping container appears damaged in any way, contact the shipping company immediately. Save all damaged packing materials to assist in proving liability for damage. Carefully inspect your exam table top while you unpack it. If any damage is noted, or if parts appear to be missing, call SSCI Customer Service at (800) 323-7366.

### Returning the Exam Table for Repairs

#### RMA Numbers:

If your exam table should require return to SSCI for repairs, discuss the problem with one of our Customer Service Representatives. Obtain an RMA number (Return Merchandise Authorization) from them before shipping the item back.

**Note:** Merchandise returned without an RMA number will not be accepted.

#### Packing and Shipment:

If you were able to keep the exam table shipping carton, repack it into the carton and staple or tape the cover securely in place. If the original shipping carton is not available, pack the table as best you can to protect it during shipment. Ship documentation with the exam table including:

- Destination
- RMA Number
- Your name, company, and address
- Your telephone number
- A description of the reason for returning the cart

### Warranty

Suburban Surgical Company, Inc. warrants the original purchaser that all equipment manufactured by Suburban Surgical Company, Inc. will be of the highest standards in material and workmanship. All equipment manufactured by Suburban Surgical Company, Inc. will be warrantied for a period of (1) year from the date of shipment from the factory.

Components and casters Suburban Surgical Company, Inc. purchases from other manufacturers will be covered by the respective manufacturer's Warranty.

Warranties will not apply if it is determined by Suburban Surgical Company, Inc. that the equipment became defective due to an accident, misuse, abuse, or alteration. Warranties do not include freight charges for replacement or repair.

## Chapter 2 – Assembly & Installation

### General

This section guides you in assembling and installing the SSCI Classic Folding Wall-Mounted Exam Table. If you should require additional assistance, please feel free to call SSCI Customer Service at 1-800-323 7366.

### Pre-assembly Preparation

On new construction, we recommend that, prior to finishing the wall to which this table will be mounted, you exactly determine where on the wall it will be located. Back up the mounting screw locations with 2" x 6"s. On an existing wall, you may want to open the wall and install these supports. These steps are not required for masonry walls.

Use lag screws for mounting into wood, or toggle bolts for mounting into block.

### Assembly Procedures

**CAUTION: Unpacking and assembling the exam table is not difficult. However, the table is heavy and we recommend that unpacking and assembly be done by at least two people. Follow the instructions carefully to prevent injury to yourself or damage to the table.**

### Tools and Supplies Required

You will need the following tools and supplies to install your wall-mounted table:


- Pen or pencil
- Tape measure
- Carpenter's Level
- Power drill with appropriate size bits
- Mounting hardware suitable for wall
- Wrenches or screwdrivers for mounting bolt/screw

### Drilling the Mounting Holes

1. Locate and drill the ten mounting holes as follows. Refer to Figure 2-2 for dimensions. Select a drill bit size suitable for your selected mounting hardware.

**Note:** The 16.948-in. height shown in Figure 2-2 is our recommendation and will result in a table height of about 42 in. You can mount the table higher or lower if it is more convenient for you.

2. Drill these ten holes into the wall.



Mounting Dimensions

Figure 2-1

### Mounting the Table to the Wall

Mount the table to the wall as follows. Refer to Figure 2-5.

**CAUTION:** Mounting the folding table to the wall is not difficult. However, the table is heavy and mounting it is very awkward for one person. We recommend that mounting be done by at least two people. Follow the instructions carefully to prevent injury to yourself or damage to the table.

1. Remove the table from the shipping carton and lay it on the floor with the table top surface down (Figure 2-2).  
**Note:** Place sheets or other soft materials on the floor to protect the table top.
2. Pull on outward on the locking handle and unfold the table (Figure 2-3).
3. Using the hardware suitable to the type of wall, mount the table to the wall with screws or bolts through the ten holes (Figure 2-4).





Figure 2-2




Figure 2-3




Figure 2-4

***Installation of your new SSCI Lateral Wall-Mounted Exam Table is now complete.***

**Chapter 3 – Folding and Unfolding the Table**

**Folding the Table**

To unfold the table and place it in the lowered position, pull up on the locking handle (Figure 3-1) and gently lower the table.




**Figure 3-1**

**CAUTION:** This table is designed to carry weights up to, but not exceeding, 200-lbs. Placing weights greater than 200-lbs on the table can damage the table or create a hazardous situation.

**Unfolding the Table**

To fold the table and place in the stored position, pull out on the locking handle (Figure 3-2) and gently lift the table up against the wall until it locks.

**Figure 3-2**



**Locking Handle**

## Chapter 4 – Accessory (Tie-down Brackets)

### Overview

A set of tie-down brackets is available for use with your exam table (**Order P/N 12805-00-AAAAAA**). The set includes four tie-down brackets and 12 Phillips-head, self-tapping screws (Figure 3-1). To install these tie-downs, follow the steps below.

Two tie-downs should be installed to each side of the table.

**Note:** The tie-downs mount to the edge or the top of the table, however, the edge is the preferred location. The locations of the tie-downs depend on the size of the animals to be examined. For smaller animals, the tie-downs should be installed close to the center of the table. For larger animals, the tie-downs should be installed farther apart, closer to the ends of the table.

### Tools Required:

- Tape measure or ruler
- Marking tool (pencil, pen, ect.)
- Punch
- Hammer
- Electric drill
- 3/32-inc. drill bit
- Phillips head screwdriver (or drive for electric drill)

**CAUTION:** Before drilling holes, tape the tie-downs where you intend to install them. Raise and lower the table to make sure that they do not interfere with the table's ability to fold up, or lie flat against the wall.

### Procedure

1. Determine where on the exam table you want the tie-downs installed.
  - Note:** Avoid placing the brackets over the existing screw heads in the of the table.
2. Hold a single tie-down in place on the edge of the table (Figure 3-2) or the top of the table (Figure 3-3), and mark the centers of the three holes on the table. The single, center hole should go toward the bottom, or toward the edge.

**CAUTION:** Before drilling holes, be very sure that the tie-downs will not interfere with the table's ability to fold up or lie flat against the wall. **DO NOT** mount a tiedown where it interferes with the table's motion or storage.

3. With a punch and hammer, tap a small dimple in the metal at each of the hole center marks. These dimples serve as guides (pilot holes) for the drill.
4. Drill the three holes with a 3/32 in. drill bit.
5. Mount the tie-down to the table with three of the screws provided in the kit.
6. Mount the remaining three tie-downs the same way.

Tie-down Bracket Kit Contents




Figure 3-1



Figure 3-2



Figure 3-3

## Chapter 5 – Cleaning

### Introduction

You will no doubt want to clean your exam table whenever it becomes dirty or saturated with waste fluids. Maintaining high standards of sanitation will be an important priority for your facility.

### Cleaning Procedures

Whenever necessary, rinse the table with clear water and dry thoroughly with clean, soft cloths. Ordinary deposits of waste and fluids can usually be removed with soap and water. Stubborn deposits may require scrubbing with “stainless steel” wool, nylon or plastic scrubbers and/or the use of commercial cleaning products. Always scrub in the direction of the “grain” of the metal. Rinse with clear water and dry thoroughly with clean, soft cloths.

Minor scale build-up and some hard water spotting may be removed by washing with vinegar, followed by a neutralizing rinse of clear water and a thorough drying with clean, soft cloths.

For heavy deposits of scale, 5% oxalic acid (use warm), 5-15% sulfamic acid, or 5-10% phosphoric acid may be used. As always, rinse with clear water and dry thoroughly with clean soft cloths.

Avoid prolonged use of chlorides (such as chlorine bleach), bromides, iodides, and thiocyanates. Never allow salty solutions to dry on the stainless steel.

**CAUTION: The warranty for this product is void if the care and cleaning instructions provided in this manual are not followed.**

© Copyright 2013 by **Suburban Surgical Co., Inc.** All rights reserved.

No part of this document may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying, recording, or by any information storage and retrieval system without written permission.

Inquiries should be addressed to **Suburban Surgical Co., Inc.**  
Wheeling, Illinois 60090